



SHC3 Antibody, FITC conjugated

Product Code	CSB-PA838779LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92529
Immunogen	Recombinant Human SHC-transforming protein 3 protein (323-468AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SHC-transforming protein 3 (Neuronal Shc) (N-Shc) (Protein Rai) (SHC-transforming protein C) (Src homology 2 domain-containing-transforming protein C3) (SH2 domain protein C3), SHC3, NSHC SHCC
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SHC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.