



UBE2E3 Antibody, Biotin conjugated

Product Code	CSB-PA839280LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969T4
Immunogen	Recombinant Human Ubiquitin-conjugating enzyme E2 E3 protein (2-128AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-conjugating enzyme E2 E3 (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme E3) (UbcH9) (Ubiquitin carrier protein E3) (Ubiquitin-conjugating enzyme E2-23 kDa) (Ubiquitin-protein ligase E3), UBE2E3, UBCE4 UBCH9
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	UBE2E3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.