



SCG3 Antibody, Biotin conjugated

Product Code	CSB-PA845166LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXD2
Immunogen	Recombinant Human Secretogranin-3 protein (205-432AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secretogranin-3 (Secretogranin III) (SgIII), SCG3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SCG3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.