



ELMO3 Antibody, Biotin conjugated

Product Code	CSB-PA846590LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BJ8
Immunogen	Recombinant Human Engulfment and cell motility protein 3 protein (464-596AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Engulfment and cell motility protein 3, ELMO3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ELMO3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.