




Pdk3 Antibody

Product Code	CSB-PA852852LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q922H2
Immunogen	Recombinant Mouse [Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 3, mitochondrial protein (131-362AA)
Raised In	Rabbit
Species Reactivity	Mouse, Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	[Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 3, mitochondrial (EC 2.7.11.2) (Pyruvate dehydrogenase kinase isoform 3), Pdk3
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Pdk3

Image



Western Blot

Positive WB detected in: Rat liver tissue

All lanes: Pdk3 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 48 kDa

Observed band size: 48 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.