



Pdk3 Antibody, Biotin conjugated

Product Code	CSB-PA852852LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q922H2
Immunogen	Recombinant Mouse [Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 3, mitochondrial protein (131-362AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	[Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 3, mitochondrial (EC 2.7.11.2) (Pyruvate dehydrogenase kinase isoform 3), Pdk3
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Pdk3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.