



ARPC1A Antibody, FITC conjugated

Product Code	CSB-PA856421LC01HU
Abbreviation	Actin-related protein 2/3 complex subunit 1A
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92747
Immunogen	Recombinant Human Actin-related protein 2/3 complex subunit 1A protein (175-317AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Probably functions as component of the Arp2/3 complex which is involved in regulation of actin polymerization and together with an activating nucleation-promoting factor (NPF) mediates the formation of branched actin networks.
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin-related protein 2/3 complex subunit 1A (SOP2-like protein), ARPC1A, SOP2L
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARPC1A