



NELL2 Antibody, HRP conjugated

Product Code	CSB-PA859098LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99435
Immunogen	Recombinant Human Protein kinase C-binding protein NELL2 protein (511-637AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein kinase C-binding protein NELL2 (NEL-like protein 2) (Nel-related protein 2), NELL2, NRP2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NELL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.