



HSD3B7 Antibody, Biotin conjugated

Product Code	CSB-PA861999LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2F3
Immunogen	Recombinant Human 3 beta-hydroxysteroid dehydrogenase type 7 protein (131-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3 beta-hydroxysteroid dehydrogenase type 7 (3 beta-hydroxysteroid dehydrogenase type VII) (3-beta-HSD VII) (3-beta-hydroxy-Delta(5)-C27 steroid oxidoreductase) (C(27) 3-beta-HSD) (EC 1.1.1.-) (Cholest-5-ene-3-beta,7-alpha-diol 3-beta-dehydrogenase) (EC 1.1.1.181), HSD3B7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSD3B7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.