



PREB Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA867205LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9HCU5 |
| Immunogen | Recombinant Human Prolactin regulatory element-binding protein (98-283AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Prolactin regulatory element-binding protein (Mammalian guanine nucleotide exchange factor mSec12), PREB, SEC12 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | PREB |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |