




TINF2 Antibody

Product Code	CSB-PA874799LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BSI4
Immunogen	Recombinant Human TERF1-interacting nuclear factor 2 protein (210-451AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TERF1-interacting nuclear factor 2 (TRF1-interacting nuclear protein 2), TINF2, TIN2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TINF2

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA874799LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.