




NDUFA4L2 Antibody

Product Code	CSB-PA878874LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRX3
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 4-like 2 protein (48-87AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 4-like 2 (NADH-ubiquinone oxidoreductase MLRQ subunit homolog) (NUOMS), NDUFA4L2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NDUFA4L2

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA878874LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.