




AK3 Antibody

Product Code	CSB-PA883438LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIJ7
Immunogen	Recombinant Human GTP:AMP phosphotransferase AK3, mitochondrial protein (142-227AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTP:AMP phosphotransferase AK3, mitochondrial (EC 2.7.4.10) (Adenylate kinase 3) (AK 3) (Adenylate kinase 3 alpha-like 1), AK3, AK3L1 AK6 AKL3L
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AK3

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA883438LA01HU at dilution of 1:100