



PRDM13 Antibody, Biotin conjugated

Product Code	CSB-PA887981LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4Q3
Immunogen	Recombinant Human PR domain zinc finger protein 13 protein (2-255AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PR domain zinc finger protein 13 (EC 2.1.1.-) (PR domain-containing protein 13), PRDM13, PFM10
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PRDM13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.