



PCDHB16 Antibody, HRP conjugated

Product Code	CSB-PA889094LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRJ7
Immunogen	Recombinant Human Protocadherin beta-16 protein (215-485AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protocadherin beta-16 (PCDH-beta-16) (Protocadherin-3X), PCDHB16, KIAA1621 PCDH3X
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PCDHB16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.