



DNAJB11 Antibody, FITC conjugated

Product Code	CSB-PA890663LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBS4
Immunogen	Recombinant Human DnaJ homolog subfamily B member 11 protein (53-238AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DnaJ homolog subfamily B member 11 (APOBEC1-binding protein 2) (ABBP-2) (DnaJ protein homolog 9) (ER-associated DNAJ) (ER-associated Hsp40 co-chaperone) (Endoplasmic reticulum DNA J domain-containing protein 3) (ER-resident protein ERdj3) (ERdj3) (ERj3p) (HEDJ) (Human DnaJ protein 9) (hDj-9) (PWP1-interacting protein 4), DNAJB11, EDJ ERJ3 HDJ9
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DNAJB11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.