



PCDHGC3 Antibody, Biotin conjugated

Product Code	CSB-PA892155LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UN70
Immunogen	Recombinant Human Protocadherin gamma-C3 protein (24-60AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protocadherin gamma-C3 (PCDH-gamma-C3) (Protocadherin-2) (Protocadherin-43) (PC-43), PCDHGC3, PCDH2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PCDHGC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.