



NDUFB9 Antibody, FITC conjugated

Product Code	CSB-PA897317LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6M9
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 9 protein (50-179AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 9 (Complex I-B22) (CI-B22) (LYR motif-containing protein 3) (NADH-ubiquinone oxidoreductase B22 subunit), NDUFB9, LYRM3 UQOR22
Immunogen Species	Homo sapiens (Human)
Research Area	Metabolism
Target Names	NDUFB9