



F11R Antibody, FITC conjugated

Product Code	CSB-PA897579LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y624
Immunogen	Recombinant Human Junctional adhesion molecule A protein (31-224AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Junctional adhesion molecule A (JAM-A) (Junctional adhesion molecule 1) (JAM-1) (Platelet F11 receptor) (Platelet adhesion molecule 1) (PAM-1) (CD antigen CD321), F11R, JAM1 JCAM
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	F11R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.