




Recombinant Human Mitochondrial import inner membrane translocase subunit Tim21 (TIMM21)

Product Code	CSB-EP871586HU
Relevance	Participates in the translocation of transit peptide-containing proteins across the mitochondrial inner membrane. Also required for assembly of mitochondrial respiratory chain complex I and complex IV as component of the MITRAC (mitochondrial translation regulation assembly intermediate of cytochrome c oxidase complex) complex. Probably shuttles between the presequence translocase and respiratory-chain assembly intermediates in a process that promotes incorporation of early nuclear-encoded subunits into these complexes.
Abbreviation	Recombinant Human TIMM21 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BVV7
Alias	TIM21-like protein, mitochondrial
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MICTFLRAVQYTEKLHRSSAKRLLLPYIVLNKACLKTEPSLRCLQYQKKTLRP RCILGVTQKTIWTQGPSRKAKEDGSKQVSVHRSQRGGTAVPTSQKVKEAGR DFTYLIVVLFGISITGGLFYTIFKELFSSSSPSKIYGRALEKCRSHPEVIGVFGESV KGYGEVTRRGRQRHVRFTFYVKDGLKHTCVKFYIEGSEPGKQGTVYAQVKEN PGSGEYDFRYIFVEIESYPRRTIIIEDNRSQDD
Research Area	Signal Transduction
Source	E.coli
Target Names	TIMM21
Expression Region	1-248aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	55.2kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.