




Recombinant Human Tropomodulin-4 (TMOD4)

Product Code	CSB-EP609953HU
Relevance	Blocks the elongation and depolymerization of the actin filaments at the pointed end. The Tmod/TM complex contributes to the formation of the short actin protofilament, which in turn defines the geometry of the membrane skeleton.
Abbreviation	Recombinant Human TMOD4 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NZQ9
Alias	Skeletal muscle tropomodulin
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MSSYQKELEKYRDIIDEILRTLSPEELEQLDCELQEMDPENMLLPAGLRQRD QTKKSPTGPLDREALLQYLEQQALEVKERDDLVPFTGEKKGKPYIQPKREIPAE EQITLEPELEEALAHATDAEMCDIAAILDMYTLMSNKQYYDALCSGEICNTEGIS SVVQPDKYKVPDEPPNPTNIEEILKRVRSDNDKEEEVNLNNIQDIPMLSELC EAMKANTYVRSFSLVATRSGDPIANAVADMLRENRLQSLNIESNFISSTGLMA VLKAVRENATLTRELVDNQRQWPGDAVEMEMATVLEQCPSIVRFGYHFTQQG PRARAAQAMTRNNELRRQQKKR
Research Area	Signal Transduction
Source	E.coli
Target Names	TMOD4
Expression Region	1-345aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	66.3kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.