




Recombinant Rat Lymphotoxin protein (Xcl1) (Active)

| | |
|----------------------------|---|
| Product Code | CSB-AP001501RA |
| Abbreviation | Recombinant Rat Ltn protein (Active) |
| Uniprot No. | P51672 |
| Form | Lyophilized powder |
| Storage Buffer | Lyophilized from a 0.2 µm filtered PBS, pH 7.4 |
| Product Type | Chemokine |
| Immunogen Species | Rattus norvegicus (Rat) |
| Biological Activity | Fully biologically active when compared to standard. The ED50 as determined by a chemotaxis bioassay using human XCR1 transfected murine BaF3 cells is less than 100 ng/ml, corresponding to a specific activity of >1.0x10 ⁴ IU/mg. |
| Purity | ≥97% as determined by SDS-PAGE. |
| Sequence | VGTEVLQESI CVSLRTQRLP VQKIKTYTIK EGAMRAVIFV TKRGLRICAD PQAKWVKTAI KTVDGRASAS KSKAETIPTQ AQRSASTAVT LTG |
| Research Area | Immunology |
| Source | E.Coli |
| Target Names | Xcl1 |
| Expression Region | 22-114aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag-Free |
| Mol. Weight | 10.0 kDa |
| Protein Length | Full Length of Mature Protein |
| PubMed ID | 9717977 |

Image



**Endotoxin**

Less than 1.0 EU/ μ g as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.