



BCAR1 Antibody, FITC conjugated

Product Code	CSB-PA002593LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56945
Immunogen	Recombinant Human Breast cancer anti-estrogen resistance protein 1 protein (68-218AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Breast cancer anti-estrogen resistance protein 1 (CRK-associated substrate) (Cas scaffolding protein family member 1) (p130cas), BCAR1, CAS CASS1 CRKAS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BCAR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.