



C1QC Antibody, HRP conjugated

Product Code	CSB-PA003641LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02747
Immunogen	Recombinant Human Complement C1q subcomponent subunit C protein (27-157AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement C1q subcomponent subunit C, C1QC, C1QG
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	C1QC