



C1QC Antibody, Biotin conjugated

Product Code	CSB-PA003641LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02747
Immunogen	Recombinant Human Complement C1q subcomponent subunit C protein (27-157AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement C1q subcomponent subunit C, C1QC, C1QG
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	C1QC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.