



CD33 Antibody, Biotin conjugated

Product Code	CSB-PA004925LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20138
Immunogen	Recombinant Human Myeloid cell surface antigen CD33 protein (18-259AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myeloid cell surface antigen CD33 (Sialic acid-binding Ig-like lectin 3) (Siglec-3) (gp67) (CD antigen CD33), CD33, SIGLEC3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD33
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.