



CD5 Antibody, FITC conjugated

Product Code	CSB-PA004942LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06127
Immunogen	Recombinant Human T-cell surface glycoprotein CD5 protein (28-234AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell surface glycoprotein CD5 (Lymphocyte antigen T1/Leu-1) (CD antigen CD5), CD5, LEU1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.