



CLDN5 Antibody, FITC conjugated

Product Code	CSB-PA005507LC01HU
Abbreviation	Claudin-5
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00501
Immunogen	Recombinant Human Claudin-5 protein (181-218AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Plays a major role in tight junction-specific obliteration of the intercellular space.
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Claudin-5 (Transmembrane protein deleted in VCFS) (TMDVCF), CLDN5, AWAL TMVCF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CLDN5