



DCLK1 Antibody, FITC conjugated

Product Code	CSB-PA006548LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15075
Immunogen	Recombinant Human Serine/threonine-protein kinase DCLK1 protein (143-381AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase DCLK1 (EC 2.7.11.1) (Doublecortin domain-containing protein 3A) (Doublecortin-like and CAM kinase-like 1) (Doublecortin-like kinase 1), DCLK1, DCAMKL1 DCDC3A KIAA0369
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DCLK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.