



EEA1 Antibody, Biotin conjugated

Product Code	CSB-PA007405LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15075
Immunogen	Recombinant Human Early endosome antigen 1 protein (734-875AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Early endosome antigen 1 (Endosome-associated protein p162) (Zinc finger FYVE domain-containing protein 2), EEA1, ZFYVE2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	EEA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.