



ERAP1 Antibody, Biotin conjugated

Product Code	CSB-PA007760LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZ08
Immunogen	Recombinant Human Endoplasmic reticulum aminopeptidase 1 protein (642-941AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoplasmic reticulum aminopeptidase 1 (EC 3.4.11.-) (ARTS-1) (Adipocyte-derived leucine aminopeptidase) (A-LAP) (Aminopeptidase PILS) (Puromycin-insensitive leucyl-specific aminopeptidase) (PILS-AP) (Type 1 tumor necrosis factor receptor shedding aminopeptidase regulator), ERAP1, APPILS ARTS1 KIAA0525
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ERAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.