



GHR Antibody, Biotin conjugated

Product Code	CSB-PA009411LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10912
Immunogen	Recombinant Human Growth hormone receptor protein (19-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth hormone receptor (GH receptor) (Somatotropin receptor) [Cleaved into: Growth hormone-binding protein (GH-binding protein) (GHBP) (Serum-binding protein)], GHR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GHR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.