



HSF1 Antibody, Biotin conjugated

Product Code	CSB-PA010791LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00613
Immunogen	Recombinant Human Heat shock factor protein 1 protein (420-529AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock factor protein 1 (HSF 1) (Heat shock transcription factor 1) (HSTF 1), HSF1, HSTF1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	HSF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.