



IL2RA Antibody, FITC conjugated

Product Code	CSB-PA011649LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01589
Immunogen	Recombinant Human Interleukin-2 receptor subunit alpha protein (22-185AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-2 receptor subunit alpha (IL-2 receptor subunit alpha) (IL-2-RA) (IL-2R subunit alpha) (IL2-RA) (TAC antigen) (p55) (CD antigen CD25), IL2RA
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL2RA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.