



INSIG1 Antibody, Biotin conjugated

Product Code	CSB-PA011745LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15503
Immunogen	Recombinant Human Insulin-induced gene 1 protein (1-87AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Insulin-induced gene 1 protein (INSIG-1), INSIG1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	INSIG1