



ITGA5 Antibody, HRP conjugated

Product Code	CSB-PA011868LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08648
Immunogen	Recombinant Human Integrin alpha-5 protein (220-472AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-5 (CD49 antigen-like family member E) (Fibronectin receptor subunit alpha) (Integrin alpha-F) (VLA-5) (CD antigen CD49e) [Cleaved into: Integrin alpha-5 heavy chain; Integrin alpha-5 light chain], ITGA5, FNRA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITGA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.