



KLK2 Antibody, Biotin conjugated

Product Code	CSB-PA012453LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20151
Immunogen	Recombinant Human Kallikrein-2 protein (25-261AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kallikrein-2 (EC 3.4.21.35) (Glandular kallikrein-1) (hGK-1) (Tissue kallikrein-2), KLK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KLK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.