




MgII Antibody

Product Code	CSB-PA013787LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8R431
Immunogen	Recombinant Rat Monoglyceride lipase protein (1-303AA)
Raised In	Rabbit
Species Reactivity	Rat, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monoglyceride lipase (MGL) (EC 3.1.1.23) (Monoacylglycerol lipase) (MAGL), MgII, Mgl2
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	MgII

Image



Western Blot

Positive WB detected in: Rat lung tissue, Mouse brain tissue

All lanes: MgII antibody at 2.8µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 34 kDa

Observed band size: 34 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.