



MgII Antibody, Biotin conjugated

Product Code	CSB-PA013787LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8R431
Immunogen	Recombinant Rat Monoglyceride lipase protein (1-303AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monoglyceride lipase (MGL) (EC 3.1.1.23) (Monoacylglycerol lipase) (MAGL), MgII, Mgl2
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	MgII
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.