




Pcsk2 Antibody

Product Code	CSB-PA017642LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28841
Immunogen	Recombinant Rat Neuroendocrine convertase 2 protein (109-637AA)
Raised In	Rabbit
Species Reactivity	Rat, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuroendocrine convertase 2 (NEC 2) (EC 3.4.21.94) (KEX2-like endoprotease 2) (Prohormone convertase 2) (Proprotein convertase 2) (PC2), Pcsk2, Nec-2 Nec2 Rpc2
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Pcsk2

Image



Western Blot

Positive WB detected in: Jurkat whole cell lysate

All lanes: Pcsk1 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 71 kDa

Observed band size: 71 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.