



# Pcsk2 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA017642LB01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P28841
<b>Immunogen</b>	Recombinant Rat Neuroendocrine convertase 2 protein (109-637AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Neuroendocrine convertase 2 (NEC 2) (EC 3.4.21.94) (KEX2-like endoprotease 2) (Prohormone convertase 2) (Proprotein convertase 2) (PC2), Pcsk2, Nec-2, Nec2 Rpc2
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	Pcsk2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.