




Pdpx Antibody

Product Code	CSB-PA017749LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60487
Immunogen	Recombinant Mouse Pyridoxal phosphate phosphatase protein (1-292AA)
Raised In	Rabbit
Species Reactivity	Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyridoxal phosphate phosphatase (PLP phosphatase) (EC 3.1.3.3) (EC 3.1.3.74) (Chronophin), Pdpx, Cin Plp Plpp
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Pdpx

Image



Western Blot

Positive WB detected in: Rat brain tissue

All lanes: Pdpx antibody at 2.8µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 32 kDa

Observed band size: 32 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.