



Pdpx Antibody, HRP conjugated

Product Code	CSB-PA017749LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60487
Immunogen	Recombinant Mouse Pyridoxal phosphate phosphatase protein (1-292AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyridoxal phosphate phosphatase (PLP phosphatase) (EC 3.1.3.3) (EC 3.1.3.74) (Chronophin), Pdpx, Cin Plp Plpp
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Pdpx
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.