





PROKR2 Antibody

Product Code	CSB-PA018749YA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFJ6
Immunogen	Recombinant Human Prokineticin receptor 2 protein (1-45AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prokineticin receptor 2 (PK-R2) (G-protein coupled receptor 73-like 1) (G-protein coupled receptor I5E) (GPR73b) (GPRg2), PROKR2, GPR73L1 PKR2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PROKR2

Image



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA018749YA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA018749YA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.