



PROKR2 Antibody, Biotin conjugated

Product Code	CSB-PA018749YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFJ6
Immunogen	Recombinant Human Prokineticin receptor 2 protein (1-45AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prokineticin receptor 2 (PK-R2) (G-protein coupled receptor 73-like 1) (G-protein coupled receptor I5E) (GPR73b) (GPRg2), PROKR2, GPR73L1 PKR2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PROKR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.