



REST Antibody, FITC conjugated

Product Code	CSB-PA019571LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13127
Immunogen	Recombinant Human RE1-silencing transcription factor protein (250-350AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RE1-silencing transcription factor (Neural-restrictive silencer factor) (X2 box repressor), REST, NRSF XBR
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	REST