



SFPQ Antibody, HRP conjugated

Product Code	CSB-PA021136LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23246
Immunogen	Recombinant Human Splicing factor, proline- and glutamine-rich protein (499-598AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Splicing factor, proline- and glutamine-rich (100 kDa DNA-pairing protein) (hPOMP100) (DNA-binding p52/p100 complex, 100 kDa subunit) (Polypyrimidine tract-binding protein-associated-splicing factor) (PSF) (PTB-associated-splicing factor), SFPQ, PSF
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SFPQ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.