



SIGMAR1 Antibody, FITC conjugated

Product Code	CSB-PA021320LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99720
Immunogen	Recombinant Human Sigma non-opioid intracellular receptor 1 protein (112-169AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sigma non-opioid intracellular receptor 1 (Aging-associated gene 8 protein) (SR31747-binding protein) (SR-BP) (Sigma 1-type opioid receptor) (SIG-1R) (Sigma1-receptor) (Sigma1R) (hSigmaR1), SIGMAR1, OPRS1 SRBP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SIGMAR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.