



SLC39A8 Antibody, HRP conjugated

Product Code	CSB-PA021637LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C0K1
Immunogen	Recombinant Human Zinc transporter ZIP8 protein (30-132AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter ZIP8 (BCG-induced integral membrane protein in monocyte clone 103 protein) (LIV-1 subfamily of ZIP zinc transporter 6) (LZT-Hs6) (Solute carrier family 39 member 8) (Zrt- and Irt-like protein 8) (ZIP-8), SLC39A8, BIGM103 ZIP8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC39A8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.