



SYNM Antibody, Biotin conjugated

Product Code	CSB-PA023019LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15061
Immunogen	Recombinant Human Synemin protein (1147-1417AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synemin (Desmuslin), SYNM, DMN KIAA0353 SYN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SYNM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.